





Installation on the sauter OFL1.0

Variant 1: Installation through the pre-drilled opening on OFL1.0






Guide the scale through the pre-drilled opening on the column for height adjustment.



Fixing of the measuring scale on the underside of the insert plate. The holes for this are already pre-drilled.

Installation on the sauter OFL1.0

Variant 2: Attaching to the holding device of the mounting plate



Mounting of the measuring scale on the underside of the inlay plate.



Screw to the holding device of the mounting plate. Required holes in OFL 1.0 are already pre-drilled. If necessary, the length of the magnetic holder must be altered.